

# Determination of 11 cannabinoids in Chocolate using \$14,990 HPLC from CTInstruments

Accurate determination of cannabinoids in chocolate is important from pricing, quality assurance, and regulatory compliance point of view. We present an easy-to-use, accurate, reliable, and affordable HPLC for measuring 11 cannabinoids in a variety of samples. This application note describes analysis of cannabis edibles.

## HPLC Features

- Reciprocating Pump
- Rheodyne 7725i Injector
- CTI HPLC Software
- UV/VIS Detector
- Temperature-controlled Column Compartment

## HPLC Specifications

<b>Flow Rate</b>	0.001 - 5mL/min
<b>Max Pressure</b>	6,300 psi
<b>Flow Accuracy</b>	±1%
<b>Flow Precision</b>	RSD <0.1%
<b>Qualitative Repeatability</b>	RSD ≤0.2% (Naphthalene/ Methanol standards)
<b>Quantitative Repeatability</b>	RSD ≤0.5% (Naphthalene/ Methanol standards)
<b>Wavelength Range</b>	180 - 680nm
<b>Spectrum Bandwidth</b>	8nm
<b>Wavelength Accuracy</b>	±1nm
<b>Wavelength Precision</b>	Below 0.1nm
<b>Noise</b>	≤0.25X10 <sup>-5</sup> AU

## HPLC Column Specifications

<b>Column Type</b>	C18, SS body
<b>Dimensions</b>	150x4.6mm
<b>Packing</b>	5µm particles
<b>Guard Column</b>	C18



**2020 COMMERCIAL  
 CANNABIS AWARDS  
 WINNER**

**Best Cannabinoids  
 Potency Testing  
 Solutions Manufacturer -  
 North America**



**cannabishplcanalyzer.com**

## Sample Information

Sample Type	Chocolate
Brand	Bhang CBD Milk Chocolate 10.00 g
Total THC per Unit	0.50 mg
Total CBD per Unit	10.00 mg



## PROCESS

### 1. Extraction

Extraction of cannabinoids from chocolate is the initial step in the analysis.

#### Extraction Parameters


Sample Weight	70 mg
Sample Preparation	grating & sonication
Extraction Solvent	methanol
Extraction Conditions	room temperature
Dilution	in acetonitrile

### 2. Injection and HPLC Analysis

After the extraction is completed, diluted extract is injected into HPLC for analysis.

#### Chromatographic Conditions


Mode	Isocratic
Temperature	30°C
Detection	UV at 220nm
Mobile Phase	Buffer:Acetonitrile
Flow Rate	1.2mL/min



### 3. Report Generation

After the analysis is completed, CTI HPLC software auto-processes the chromatogram, followed by export to custom lab report generation program in MS Excel (highly customizable and automated report generation for ease of use).

#### CANNABINOID PROFILE



Compound		Result (% w/w)	mg/gram of sample
THCV	Tetrahydrocannabivarin	NR	NR
Δ8-THC	(-)-Δ8-THC	NR	NR
Δ9-THC	(-)-Δ9-THC	<0.01	<0.05
Δ9-THCA-A	(-)-trans-Δ9-THC acid A	NR	NR
CBD	Cannabidiol	0.11	1.07
CBDA	Cannabidiolic acid	NR	NR
CBDV	Cannabidivarin	NR	NR
CBG	Cannabigerol	NR	NR
CBGA	Cannabigerolic acid	NR	NR
CBN	Cannabinol	NR	NR
CBC	(+/-) Cannabichromene	NR	NR
<b>Total Cannabinoids *</b>		<b>0.11</b>	<b>1.07</b>
Total Potential THC		<0.01	0.04
Total Potential CBD		0.11	1.07
Total Potential CBG		<0.01	<0.05

#### Results

	Manufacturer's Values	Measured Values
<b>Total THC per Unit</b>	0.50 mg	0.40 mg
<b>Total CBD per Unit</b>	10.00 mg	10.70 mg

#### Lower Limit of Quantification (LLOQ)

The lower limit of quantification (LLOQ) is the lowest amount of a cannabinoid in a sample that can be quantitatively determined with suitable precision and accuracy using the corresponding method and dilution rates. All values below this threshold are reported as NR - None Reported.

Compound		LLOQ (% w/w)
THCV	Tetrahydrocannabivarin	0.01
Δ8-THC	(-)-Δ8-THC	0.01
Δ9-THC	(-)-Δ9-THC	0.01
Δ9-THCA-A	(-)-trans-Δ9-THC acid A	0.01
CBD	Cannabidiol	0.01
CBDA	Cannabidiolic acid	0.01
CBDV	Cannabidivarin	0.01
CBG	Cannabigerol	0.01
CBGA	Cannabigerolic acid	0.01
CBN	Cannabinol	0.01
CBC	(+/-) Cannabichromene	0.01

#### Instrument Calibration & Quality Control

Date of Quality Control	Standard	Standard Concentration (ug/mL)	Measured Concentration (ug/mL)	Delta (%)	PASS/FAIL	Notes
28-Jan-21	Benzoic acid	1002.9	1013.0	1.0%	PASS	
28-Jan-21	CBD	100.5	103.4	2.9%	PASS	